PRODUCT **SPECIFICATIONS**



WATER FILTER SYSTEMS

th 211

FILTERED.CLEANER.BETTER.

ENVIRONMENTALLY PREFERRED:

Body Glove's water filtration systems are designed with the environment in mind by: 1) reducing the number of plastic water bottles discarded into landfills, 2) using recyclable materials in packaging and 3) providing an easy method to recycle the filtration cartridges.

RESPONSIBLE LIVING: Body Glove is a leader in the recycling revolution. Purchasing a Body Glove filter allows raw materials from the filter to be recycled, reducing landfill waste and saving resources.

Call toll free: 1-888-62-REUSE for a recycle kit with a pre-paid UPS label.

FLOW RATE - 0.75 GPM

Capacity - 1,000 Gallons







MICRO-FILTER TECHNOLOGY

Body Glove's filter cartridges use a pre-filter wrap and a patented carbon block that, in combination, help trap sediment particles, chlorine taste and odor and other contaminants. This leading-edge, engineered media helps deliver cleaner, more sediment-free water resulting in improved drinking, cleaning and cooking water.

Filter system shown includes head and cartridge. Cartridges also available separately. An accessory filter faucet is required and sold separately.



www.bodyglove.waterinc.com

NSF is a trademark of NSF International

KDF® is a registered trademark of KDF® Fluid Treatment, Inc.

WQA Certified for Sustainability is a trademark of the Water Quality Association

Body Glove, The Hand Design and Protect The Core are registered trademarks of Body Glove International, LLC.

Distributed by Water, Inc. 800-322-WATER(9283) · waterinc.com











REDUCES:

Chlorine Taste & Odor **Chloramines**

Cloudiness, Sediment, Rust**

Giardia Lamblia Cysts

Cryptosporidium Parvum Cysts

Molds and Algae** Lime Scale*

Asbestos Fibers

Oxidized Iron**

Oxidized Manganese**

Oxidized Sulfides**

0.5 Micron Particulate Reduction

Total Trihalomethanes (TTHMs)

Volatile Organic Chemicals (VOCs) Methyl Tertiary-Butyl Ether (MTBE)

0.2 Micron Reduction: Phenytoin Naproxen Estrone Ibupropen Bisphenol A

Nonyl phenol





System Tested and Certified by NSF International against NSF/ANSI Standards 42, 53 and 401 in model BG-12000 for the reduction of:

Standard 42 - Aesthetic Effects

- · Chemical Reduction
 - · Chlorine Taste & Odor
- · Chloramines
- Mechanical Filtration
- Particulate Reduction- Class I

Standard 53 - Health Effects

- · Chemical Reduction
- Lead
- · VOC
- Mechanical Filtration
- Cvsts
- Asbestos

Standard 401 - Pharmaceuticals: Emerging Compounds / incidental contaminants

- Phenytoin Naproxen
- IbupropenBisphenol A
- Nonyl phenol (pending)



3004 GLOVE WATER FILTER SYSTEMS

BG-12000

FILTERED.CLEANER.BETTER.




enhance your filtered water with **INSTANT HOT WATER**

enhance your filtered water with

BG-12000 system connects to an accessory faucet (not included) choose from many styles & finishes (EverHot Modern Series Hot & Cold faucet model shown)

> O Body Glove Water Filter Systems install easily under your sink

- · Easy Installation-instructions included.
- · Micro-filtration technology and compact design offers multiple install locations saving under-counter space, and helps provide convenient and quick filter changes.
- Engineered carbon block technology conserves water for maximum filter life.
- Internal shut-off valve means less hassle and no guess work when replacing filters.
- Keep track of water useage with a digital flow meter (included with BG12000 only)



- Contaminants and substances listed for reduction by this filtration system are not necessarily present in your water.
- WARNING: Do not use with water that is microbiologically unsafe or from an unknown source or quality without sufficient disinfection before and after the filtration system. Systems certified for cysts reduction may be used on disinfected water that may contain filterable cysts. FPA
 Establishment # 070595-CT-001. Conforms to NSF/ANSI 53 for VOC reduction. See performance data sheet for individual contaminants and reduction performance. Well water applications should have a comprehensive well analysis confirming properly disinfected and potable water prior to installation of the Body Glove system.
- Body Glove filtration capacity and performance is dependent upon incoming water conditions. The filtration system must be maintained according to the manufacturer's instructions in addition to the replacement of filtration cartridges. Replace the filter annually or sooner if the flow rate becomes slow. Capacity can vary based on local water conditions. System does not include a faucet, which must be purchased separately.
- Conforms to the NSF/ANSI Standard 42 and 53 and 401 for the reduction of the listed substances. For conditions of use, health claims certified by the California Department of Public Health, and replacement parts, see product data sheet. California Department of Public Health Certification Number: 08-1939

INTEGRATED MEMBRANE PRE-ACTIVATED CARBON TECHNOLOGY

Body Glove incorporates a patented Integrated Membrane Pre-Activated Carbon (IMPACT) Technology. This filtration technology blends the high-efficiency of carbon media with the durability of a micro-porous membrane filter into one single pharmaceutical-grade cartridge. This unique process dramatically reduces pressure drop, and increases flow rates. For the consumer, this translates into longer cartridge life while eliminating the need for additional pretreatment.

QUADSTAGE FILTRATION

The Body Glove BG-12000 Filtration System incorporates an Innovative filtering design that combines four filtering processes into one recyclable polymer filter housing. Water is first treated by KDF-55 media for bacteriostatic and scale control. Next, a Micro-Filter-Wrap more evenly traps sediment particles while maintaining overall flow rates. Water is then filtered by a patented 0.2 micron pharmaceutical-grade membrane. Lastly, the filtered water passes through a patented, high-efficiency carbon block to reduce unpleasant taste and odor concerns. The Quad-Stage Filtration System helps provide a more comprehensive drinking water filtration process.

Manufactured by 3M in the USA with globally sourced components.



Distributed by Water, Inc. 1044 East Del Amo Blvd. Carson, CA 90746

Ph: 800.322.WATER (9283) Fax: 310.885.4385 waterinc.com

BG-12000: W-603 Rev. 3-15